

CRISPR/Cas9 genome-wide sgRNA libraries for loss-of-function and gain-of-function genetic screens


Donato Tedesco¹, Mikhail Makhanov¹, Debbie Deng¹, Sylvain Baron¹, Paul Diehl¹, Dmitry Suchkov¹, Costa Frangou¹, and Alex Chenchik¹.
¹Cellecta, Inc., Mountain View, CA USA

Abstract

Genome-wide loss-of-function screening is a fundamental method to identify genes responsible for driving biological processes. Complex pooled lentiviral-based libraries expressing large numbers of genetic disruptors, such as shRNAs (RNAi) or sgRNAs (CRISPR), make large-scale cell screening practical. Although loss-of-function shRNA and sgRNA pooled library screens are similar in concept, the gene loss of function is achieved by different mechanisms (mRNA degradation with RNAi, full gene disruption with CRISPR-KO, transcriptional inhibition with CRISPRi), so some divergences are expected and indeed observed when comparing results obtained using one method versus the other. Furthermore, contrary to RNAi, CRISPR technology can be modified to activate gene expression (CRISPRa), thus enabling the use of genome-wide gain-of-function screening in gene function studies. Here we present and compare results from genetic screens in human cancer cells using CRISPR-KO, CRISPRi, CRISPRa and RNAi technologies. To investigate differences in the RNAi and CRISPR HTGS technologies, we performed parallel viability screens using sgRNA and shRNA libraries targeting the same set of ~6,300 genes with the same number of effectors for each gene. The screens were aimed at the identification of essential genes in a pair of primary syngenic PDX-derived CML cell lines. The screen results show significant but not complete overlap in the genes identified by each technology in each cell line, indicating that both approaches are effective at identifying the majority of essential genes in a cell system. However, deeper analysis revealed that a small number of genes could only be identified by either CRISPR or RNAi technology. Importantly, by combining data from two screening methodologies, a higher number of viability genes and pathways could be identified and subsequently validated by independent cell-based assays with a high confirmation rate. To assess the performance of CRISPRa and CRISPRi HTGS technologies, multiple parallel screens were performed in MDA-MB-231 cells with different sgRNA design libraries targeting the same set of ~2000 genes. Gene target identification rates achieved by the different libraries were then used to identify the best sgRNA designs for CRISPRa and CRISPRi applications.

sgRNA Scaffold Improvement

Natural CRISPR Target Complex in bacteria



Standard CRISPR single-guide RNA (sgRNA)




Improved CRISPR single-guide RNA (sgRNA)



Modified sgRNA structure increases knock-out effectiveness
• Improves robustness of screen by increasing signal (sgRNA depletion or enrichment) over random genetic variation (noise).

A-T inversion removes transcriptions stop signal of TTTT
Extension of first loop to increase Cas9 binding
Combined modifications produce Helix Extended A-T inversion ("HEAT") sgRNA scaffold

Library Construction and Screening




sgRNA Structure Improvement

4X CRISPR LIBRARY:

- ~100 genes, 8 sgRNAs/gene (including known essential genes)
- ~100 non-targeting sgRNAs (negative controls)
- ~100 intron-targeting sgRNAs (negative controls)

Each sgRNA repeated 4X:

- sgRNA-wt structure
- sgRNA-AT structure
- sgRNA-HE structure
- sgRNA-HEAT structure



sgRNA structure optimization screen






Figure 1. To identify the best performing sgRNA design, we built a sgRNA design optimization library of different handle structures, all targeting the same ~900 human genomic sites. Targeted sites included essential genes, non-essential genes, introns, and non-human control sequences. In order to identify the best performing sgRNA design, the optimization library was transduced into CAS9 expressing cell lines validated with CRISPRTest Cas9 Assay kit [2]. Samples were then collected at multiple time points post transduction. On the left, the waterfall plots display the depletion fold of positive and negative controls within the library for each design.

sgRNA structure optimization screen, validation



CRISPR vs RNAi Screens



Analysis of Controls




Figure 3. Parallel CRISPR/RNAi screens in CML-R cell line: Analysis of Controls. (a) Changes in representation (depletion) of negative and positive (known cytotoxic) controls in the shRNA and gRNA libraries was monitored over a 4 week period. (b) The standardized mean depletion difference between positive and negative controls was calculated, as a means to compare RNAi vs. CRISPR screening platform performance.

shRNA / gRNA Hit Distribution




Figure 4. Measured distribution of gRNA and shRNA hits within target genes in CML-R and CML-P cell lines. In both cell lines, about 35% of target genes in the shRNA library were identified by a single shRNA hit (for a total of ~2,300 shRNAs in the shRNA library) compared to only about 5% genes (~320 gRNAs) in the gRNA library. Single shRNA/gRNA hits are most likely off-target hits. Target genes identified by multiple hits (4 or more) were found at comparable frequencies in CRISPR and RNAi screens. Multiple shRNA/gRNA hits are most likely on-target hits.

RNAi / CRISPR Gene Hit Overlap




Figure 5. Gene hit analysis: Overlap and Reproducibility. RNAi/CRISPR overlap: the number of RNAi and CRISPR gene hits are displayed for the CML-R and CML-P cell lines. The overlap between RNAi and CRISPR hits is reported in the grey bar.

RNAi / CRISPR Unique Gene Hits




Figure 6. The number and lists of gene hits exclusively identified by either CRISPR or RNAi technology (no signal with other technology) are displayed.

Pathway Analysis

| Pathway | # genes | Z Score | Adjusted p | Pathway | # genes | Z Score | Adjusted p |
|--|---------|---------|------------|--|---------|---------|------------|
| Cytosolic tRNA aminoacylation | 16 | 14.2 | 0.005 | purine nucleotides de novo biosynthesis II | 5 | 7.4 | 0.005 |
| Activation of APC/C and APC/C mediated degradation of mitotic proteins | 12 | 12.2 | 0.005 | Cap-dependent Translation Initiation | 2 | 7.3 | 0.005 |
| Regulation of APC/C activators between G2/S and early anaphase | 12 | 11.8 | 0.005 | Transferrin endocytosis and recycling | 4 | 6.6 | 0.005 |
| CD1 association with the CDC6/ORA1 origin complex | 10 | 11.1 | 0.005 | Pyrimidine biosynthesis | 2 | 6.4 | 0.005 |
| VPI-mediated degradation of APOBEC3G | 9 | 11.0 | 0.005 | Polo-like kinase mediated events | 2 | 6.4 | 0.005 |
| SCF-beta-NCP mediated degradation of Em1 | 9 | 11.0 | 0.005 | Formation of RNA Pol I elongation complex | 6 | 6.2 | 0.005 |
| mRNA Splicing - Major Pathway | 7 | 10.9 | 0.005 | Signaling by Insulin receptor | 7 | 5.9 | 0.005 |
| Cell Cycle Checkpoints | 15 | 10.7 | 0.005 | Nucleotide Excision Repair | 5 | 5.8 | 0.005 |
| Synthesis of DNA | 10 | 10.7 | 0.005 | CFR transcription factor network | 7 | 4.6 | 0.005 |
| Stabilization of p53 | 10 | 10.7 | 0.005 | superpathway of geranylgeranylphosphate biosynthesis I (via mevalonate) | 3 | 4.4 | 0.005 |
| Ork1 removal from chromatin | 10 | 10.7 | 0.005 | pyrimidine deoxyribonucleotides de novo biosynthesis | 3 | 4.4 | 0.005 |
| Autodegradation of the Ubiquitin ligase COP1 | 10 | 10.7 | 0.005 | Insulin synthesis and processing | 3 | 4.4 | 0.005 |
| Vpu mediated degradation of CD4 | 9 | 10.5 | 0.005 | ATM pathway | 18 | 4.3 | 0.005 |
| Destabilization of mRNA by AUF1 (hnRNP D0) | 9 | 10.5 | 0.005 | Validated targets of C-MYC transcriptional activation | 7 | 3.8 | 0.005 |
| Cytin A/C6/3 associated events at S phase entry | 11 | 10.4 | 0.005 | Iron uptake and transport | 4 | 4.4 | 0.009 |
| Ubiquitin-Mediated Degradation of Phosphorylated Cdc25A | 10 | 10.3 | 0.005 | Mitochondrial tRNA aminoacylation | 3 | 3.2 | 0.009 |
| Signaling by Wnt | 9 | 10.1 | 0.005 | methionine degradation I (to homocysteine) | 1 | 4.1 | 0.012 |
| SCF(Skp2)-mediated degradation of p27(p21) | 9 | 10.1 | 0.005 | Unfolded Protein Response | 5 | 2.3 | 0.012 |
| Regulation of activated PAK-2/3/4 by proteasome mediated degradation | 9 | 10.1 | 0.005 | Regulation of nuclear beta catenin signaling and target gene transcription | 5 | 2.3 | 0.012 |
| Autodegradation of beta-catenin by the destruction complex | 9 | 10.1 | 0.005 | Golgi, External Pericentriolar Stack Reorganization | 2 | 3.4 | 0.016 |
| Ubiquitin-dependent degradation of Cytin D1 | 10 | 10.1 | 0.005 | Recruitment of mitotic centrosome proteins and complexes | 2 | 3.2 | 0.016 |
| Signaling by Wnt | 13 | 9.9 | 0.005 | Centrosome maturation | 2 | 3.2 | 0.016 |
| Regulation of ornithine decarboxylase (ODC) | 9 | 9.7 | 0.005 | PI3K Cascade | 1 | 3.4 | 0.018 |
| ER-Phagosome pathway | 9 | 9.7 | 0.005 | Insulin signaling pathway | 3 | 2.8 | 0.018 |
| Mitotic G1-G2/S phases | 15 | 9.7 | 0.005 | GPCR downstream signaling | 1 | 3.4 | 0.021 |
| RNA Polymerase II Transcription | 11 | 8.6 | 0.005 | p73 transcription factor network | 4 | 2.5 | 0.021 |
| mRNA Processing | 10 | 8.5 | 0.005 | Insulin receptor signaling cascade | 2 | 2.5 | 0.021 |
| Regulation of Apoptosis | 9 | 8.4 | 0.005 | cysteine biosynthesis III (mammalia) | 1 | 3.4 | 0.024 |
| Pausing and recovery of elongation | 5 | 8.3 | 0.005 | Activation of ATR in response to replication stress | 1 | 3.1 | 0.020 |

Figure 7. Data from CRISPR and RNAi screens was combined for pathway analysis of gene hits. Identified pathways were filtered for redundancy and sorted by significance base on Z-scores and adjusted p-values (0.005<p<0.05). A subset of lower ranking pathways was then chosen for confirmation studies.

Confirmation Studies

| Pathway | Tested genes | Protein | Inhibitory Drug | Drug Toxicity |
|---|--------------|---|-----------------|---------------|
| Transferrin endocytosis and recycling / iron uptake and transport | ATRB0B | V-ATPase | Balfemycin A1 | + |
| | ATRB0C | | Balfemycin B1 | + |
| | ATRB0E | | Canacynanin A | + |
| Signaling by Insulin receptor | ATRB0D | V-ATPase | Balfemycin A1 | + |
| | ATRB0E | | Balfemycin B1 | + |
| Superpathway of geranylgeranylphosphate biosynthesis I | FDP5 | Farnesyl Diphosphate Synthase | Zileuton acid | +/+ |
| | HMGCS1 | 3-Hydroxy-3-Methylglutaryl-CoA Synthase 1 | Hymegimmon | +/+ |
| | ATRB0E | Adenosylhomocysteinase 2 | Neplirocin | + |
| Methionine degradation I (to homocysteine) | AHCY1 | S-adenosylhomocysteinase 2 | Neplirocin | + |
| | AHCY1 | Adenosylhomocysteinase 2 | Neplirocin | + |
| Cysteine biosynthesis III (mammalia) | MAT2A | Methionine Adenosyltransferase II, Alpha | FIDAS-5 | + |
| | MAT2A | Methionine Adenosyltransferase II, Alpha | FIDAS-5 | + |

Figure 8. Confirmation studies on a subset of pathways identified by combined CRISPR/RNAi data analysis were performed. The cytotoxic effect of inhibitory compounds specific to the protein products of the identified target genes was assayed in CML-R cells.

CRISPRa and CRISPRi systems





Figure 9. Cells expressing dCas9-KRAB repressor and CMV-GFP were transduced with targeting (sg_GFP) or non-targeting (sg_control) sgRNAs, expressed from either wt (U6wt) or tet-inducible (U6tet) U6 promoter. CMV repression was assayed by loss of GFP fluorescence intensity.

Figure 10. Cells expressing dCas9-VPH transactivator were transduced with sgRNAs targeting the promoters of the indicated genes. Promoter transactivation was measured by qRT-PCR.

Figure 11. Parental cells and cells expressing dCas9-VPH transactivator (both expressing GFP and RFP from CMV and UBC promoters, respectively) were transduced with UBC-targeting sgRNA. UBC promoter transactivation was assayed by fluorescence intensity.

CRISPRa/i screens



Conclusions

In the presented studies, the CRISPR/Cas9 screening platform showed comparable or better performance than RNAi.

- Modifications in the sgRNA structure enhance performance in CRISPR genetic screens
- Fewer off-target hits are seen in CRISPR vs RNAi screens
- More "pathway coverage" is seen in CRISPR vs RNAi screens
- In some instances, RNAi still performs better than CRISPR
- Combined CRISPR/RNAi screens outperform replicate CRISPR-only and RNAi-only screens
- CRISPRi-specific sgRNA design improves performance of loss-of-function CRISPR screens
- CRISPRa-specific sgRNA design improves performance of gain-of-function CRISPRa screens