



**CELLECTA**

**USER MANUAL**

# **LINE-1 Mobile Element NGS Assay**

v1 — Last update: 2018/07/04

Cellecta, Inc.

# Table of Contents

- LINE-1 Mobile Element NGS Assay ..... 1**
- LINE-1 Kit Components..... 2**
- Additional Necessary Materials ..... 4**
- Protocol Overview ..... 5**
  - Step 1. DNA Shearing ..... 7
  - Step 2. Primer Extension..... 8
  - Step 3. T-Adapter Ligation ..... 9
  - Step 4. 1st Round of PCR ..... 10
  - Step 5. 2nd Round of Nested PCR ..... 11
  - Step 6. 3rd Round PCR with Indexed Primers ..... 12
  - Step 7. Probe Preparation for Illumina NGS ..... 14
  - Step 8. Illumina NGS and Data Analysis ..... 16
- Appendices ..... 18**
  - Appendix A: Library Validation for Illumina NGS by TA Cloning (Optional Procedure) ..... 19
  - Appendix B: Structure of the 3' L1H Amplicon and Primer Locations ..... 21
  - Appendix C: Index Sequences ..... 22
- Technical Support ..... 23**
- Contact Us ..... 24**
- PDF Download Link ..... 25**

# LINE-1 Mobile Element NGS Assay

---

## Background

Two-thirds of the human genome is made up of mobile elements comprising *Long Interspersed Nuclear Elements* (LINEs, such as LINE-1). Members of the L1H subfamily of LINE-1 are active transposons that induce genetic rearrangements throughout the genome. These transpositions contribute to genetic diversity in tumor cell populations and have also been implicated in schizophrenia, Alzheimer's, and other neurodegenerative diseases. Cellecta's LINE-1 Mobile Element NGS Assay provides a convenient, efficient, and targeted NGS-based approach to identify *de novo* polymorphic insertions of LINE-1 in the human genome. The assay uses specially-designed primers to selectively amplify only the transposable elements of the youngest L1H family of LINE-1 that are still active in the genome, as opposed to the ancient elements fixed in the population. The amplified fragments are then characterized by NGS to assess the prevalence and diversity in the genome. This approach enables detection and identification of all the 3' L1H variations for a sample in a single reaction, enabling rapid analysis of multiple samples.

## Intended Use:

The LINE-1 Mobile Element NGS Assay kit is designed to profile L1H insertion sites in human DNA. This kit is for research use only.

**Please read the entire user manual before proceeding with your experiment.**

---

**PLEASE NOTE:** The purchase of all Cellecta Products is covered by Cellecta's standard Terms and Conditions of Sale as described on [Cellecta's website](#), and selected Products containing particular technology or having certain features are also subject to restrictions on use as outlined in the [Label License section](#) of our website. Please review these Terms and Label License Restrictions before opening and using your Product and, if you are not able to abide by the restriction, contact Cellecta to return the item to Cellecta for a full refund.

---



[Download the PDF Version of the LINE-1 Mobile Element NGS Assay User Manual here](#)

# LINE-1 Kit Components

---

## Components

| Description               | Cap Color  | No. of Vials | Concentration | Volume   |
|---------------------------|------------|--------------|---------------|----------|
| Taq DNA Polymerase        | pink       | 1            | 100X          | 15 µl    |
| Taq DNA Polymerase Buffer | pink       | 1            | 10X           | 150 µl   |
| GSP1* Primer              | yellow     | 1            | 10 µM         | 55 µl    |
| dNTP Mix                  | white      | 1            | 50X           | 105 µl   |
| PCR-Grade Water           | white      | 2            | NA            | 3,600 µl |
| DNA Ligase                | purple     | 1            | 100X          | 15 µl    |
| DNA Ligase Buffer         | purple     | 1            | 10X           | 130 µl   |
| T-Adapter                 | yellow     | 1            | 10 µM         | 55 µl    |
| HF DNA Polymerase         | light blue | 1            | 100X          | 40 µl    |
| HF DNA Polymerase Buffer  | light blue | 1            | 5X            | 780 µl   |
| Adapter Primer            | yellow     | 1            | 10 µM         | 55 µl    |
| GSP2 Anchor               | yellow     | 1            | 10 µM         | 30 µl    |
| P5 Primer                 | yellow     | 1            | 10 µM         | 30 µl    |
| Read 1 Seq Primer         | blue       | 1            | 100 µM        | 30 µl    |
| Index 1 Seq Primer        | blue       | 1            | 100 µM        | 30 µl    |
| Read 2 Seq Primer         | blue       | 1            | 100 µM        | 30 µl    |
| P7 Indexed Primer 1       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 2       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 3       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 4       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 5       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 6       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 7       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 8       | white      | 1            | 10 µM         | 5 µl     |
| P7 Indexed Primer 9       | white      | 1            | 10 µM         | 5 µl     |

|                      |       |   |            |           |
|----------------------|-------|---|------------|-----------|
| P7 Indexed Primer 10 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 11 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 12 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 13 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 14 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 15 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 16 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 17 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 18 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 19 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 20 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 21 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 22 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 23 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |
| P7 Indexed Primer 24 | white | 1 | 10 $\mu$ M | 5 $\mu$ l |

\*GSP = Gene-Specific Primer

**Note:** Primer sequences can be accessed on the Cellecta website.

### **Storage Conditions:**

Store kits and all components at -20°C until ready for use.

# Additional Necessary Materials

## Additional Reagents and Apparatus:

| Name   | Recommended Manufacturer | Catalog #   | Description             |
|--|--------------------------|-------------|-------------------------|
| Paired-End Cluster Generation Kit for HiSeq or NextSeq Instruments | Illumina                 | PE-402-4002 | NGS of libraries        |
| QIAquick Gel Extraction Kit  | QIAGEN                   | 28704       | DNA extraction from gel |
| QIAquick PCR Purification Kit                                      | QIAGEN                   | 28701       | PCR purification        |
| Semi-skirted 96-well PCR plates                                    | Eppendorf                | 951020303   | PCR amplification       |
| Clear Adhesive Film for 96-well plates                             | Bio-Rad                  | MSB1001     | PCR Plate seal          |

## Required Instrumentation:

The protocols were optimized using the specific instrument specified. Some modifications and optimization to the protocol may be necessary if using different models or instruments from other manufacturers.


| Name                                   | Manufacturer       | Model            |
|--|--------------------|------------------|
| Sonicator                              | Covaris            | E-series         |
| 96-Well Thermal Cycler (or equivalent) | Applied BioSystems | Veriti® 96-well  |
| Next-Gen Sequencing (NGS) Platform     | Illumina           | HiSeq®, NextSeq® |

**Other than the specific reagents and instruments above, the protocol assumes the user has access to standard materials (e.g., polypropylene tubes, pipette tips), equipment (table top centrifuges, pipettes, scales), and common reagents (e.g., TE buffer, ethanol) and buffers used in a typical life science laboratory.**

# Protocol Overview

---

The LINE-1 Mobile Element NGS Assay uses primers specifically designed to amplify the 3' L1H mobile genetic elements in randomly-sheared DNA (Step 1) to detect novel L1H genomic insertions. An L1H 3'-UTR-specific primer (GSP1) is used for a single round of extension into the chromosome (Step 2), followed by ligation of a T-Adapter to generate the template for PCR amplification (Step 3). Two rounds of amplification follow synthesis of the 3' L1H templates (Steps 4 & 5). The first round uses GSP1 and the Adapter Primer. The second round of PCR uses a nested GSP2 Anchor Primer with the Adapter Primer. A third round of amplification with primers containing the Illumina P5 and P7 adapter sequences (Step 6) prepares the samples for sequencing on the Illumina HiSeq or NextSeq instrument (Step 7).





# Step 1. DNA Shearing

---

Genomic DNA from the biological samples first needs to be sheared by sonication into small enough fragments for subsequent denaturation, primer extension, and addition of the T-Adapters. Random fragmentation of DNA by sonication does not depend on specific DNA sequences, so that the unique ends at random DNA breaks allow for detection of NGS reads coming from independent ligations to the T-Adapter which ensures reliable detection of L1H integration sites. For this reason, we do not use restriction enzymes or other sequence-dependent nucleases techniques.

**NOTE:** The protocol below was developed using the Covaris E-series sonicator. Other sonicators or mechanical techniques (e.g., shearing through a syringe) may be used to generate randomly-sheared DNA in the size range of 750-1200 bp. DNA shearing conditions need to be optimized for the specific instrument used to generate fragments of 750-1200 bp. Regardless of the instrument used for this step, fragmentation must be random. Do not use nucleases or other enzymatic approaches for the reasons noted above.

For each DNA sample to be sequenced, prepare one library.

1. Dilute 1 µg of high-quality genomic DNA (OD 260/280 ratio within a 1.8 to 2.0 range) with 1X Low TE Buffer (10:0.1) in a 1.5-ml LoBind (or similar) tube to a total volume of 50 µl.
2. Randomly fragment DNA to a size range of 750-1200 bp. For the Covaris E-series or S-series instruments, follow the manufacturer's instructions. An example of parameter settings from the Covaris E instrument are the following:

| Parameters       | Settings  |
|------------------|-----------|
| Duty factor      | 0.2       |
| Cycles per Burst | 200       |
| Treatment Time   | 6 minutes |
| Bath Temperature | 4 – 8° C  |

3. Slowly transfer sheared DNA into a PCR tube or well of the 96-well plate.

## Step 2. Primer Extension

This step primes and elongates the 3' regions of L1 elements (3' L1) using the GSP1 primer which sits ~100 bp from the start of the L1 element (Figure 2A). The Taq DNA Polymerase needs to be used in this reaction since it adds a 3'-adenine overhang to the ends of the extended double-stranded fragments. Do not use the HF DNA Polymerase, as it does not have terminal transferase activity.

1. Prepare the following reaction mix on ice in a 50  $\mu$ l thin-wall PCR tube:

| Component                               | Volume                      |
|---|-----------------------------|
| Sheared DNA sample from Step 1 (200 ng) | 10 $\mu$ l                  |
| Taq DNA Polymerase Buffer (10X)         | 5 $\mu$ l                   |
| dNTP Mix (50X)                          | 1 $\mu$ l                   |
| GSP1 Primer                             | 1 $\mu$ l                   |
| PCR-Grade Water                         | 32.5 $\mu$ l                |
| Taq DNA Polymerase (100X)               | 0.5 $\mu$ l                 |
| <b>Total</b>                            | <b>50 <math>\mu</math>l</b> |

2. Mix gently, and centrifuge briefly to collect droplets.
3. Perform primer extension using the following program:

| Temperature | Time      |
|-------------|-----------|
| 95°C        | 1 minute  |
| 62°C        | 5 minutes |
| 72°C        | 5 minutes |
| 4°C         | $\infty$  |

4. Purify the extended 3' L1 product using the QIAquick PCR Purification Kit. Elute DNA in 30  $\mu$ l of QIAGEN EB buffer.

## Step 3. T-Adapter Ligation

---

This step ligates a phosphorylated T-Adapter with a 3'-T overhang to the ends of the DNA fragments from primer extension 3' L1H products that have 3'-A overhangs generated at Step 2. The reaction adds a specific adapter sequence to the 3' ends of the genomic fragment.

1. Prepare the reaction mix on ice in the following order:

| Component                            | Volume                      |
|--------------------------------------|-----------------------------|
| Extended 3' L1H products from Step 2 | 30 $\mu$ l                  |
| T4 DNA Ligase Buffer (10X)           | 5 $\mu$ l                   |
| T-Adapter                            | 2 $\mu$ l                   |
| PCR-Grade Water                      | 12.5 $\mu$ l                |
| T4 DNA Ligase (100X)                 | 0.5 $\mu$ l                 |
| <b>Total</b>                         | <b>50 <math>\mu</math>l</b> |

2. Incubate the reaction mix for 30 minutes at 37°C.
3. Follow the instructions in the QIAquick PCR Purification Kit to purify the reaction on one QIAquick column, eluting in 30  $\mu$ l of QIAGEN EB buffer.

## Step 4. 1st Round of PCR

The First Round of PCR selectively amplifies DNA fragments that have the T-Adapter on one end and the 3' L1-specific sequences on the other. These fragments are amplified by PCR using the HF (High Fidelity) Polymerase for 15 cycles in 50  $\mu$ l using the GSP1 primer and T-Adapter primer.

1. Use half (15  $\mu$ l) of the T-Adapter-ligated 3' L1H DNA sample from Step 3. Save the other half in case you may need to repeat the experiment due to potential problems during the next steps of the protocol. Prepare the following reaction mix in a 50  $\mu$ l thin-wall PCR tube:

| Component                         | Volume                      |
|-----------------------------------|-----------------------------|
| T-Adapter-ligated DNA from Step 3 | 15 $\mu$ l                  |
| HF DNA Polymerase Buffer (5X)     | 10 $\mu$ l                  |
| dNTP Mix (50X)                    | 1 $\mu$ l                   |
| GSP1 Primer                       | 1 $\mu$ l                   |
| Adapter Primer                    | 1 $\mu$ l                   |
| PCR-Grade Water                   | 21.5 $\mu$ l                |
| HF DNA Polymerase (100X)          | 0.5 $\mu$ l                 |
| <b>Total</b>                      | <b>50 <math>\mu</math>l</b> |

2. Mix gently and centrifuge briefly to collect droplets.
3. Amplify adapter-ligated 3' L1H DNA products using the following PCR program:

|                  |                  |
|------------------|------------------|
| 98°C, 60 seconds | <b>1 cycle</b>   |
| 98°C, 15 seconds |                  |
| 62°C, 15 seconds | <b>15 cycles</b> |
| 72°C, 30 seconds |                  |
| 72°C, 60 seconds | <b>1 cycle</b>   |
| 4°C, $\infty$    | <b>1 cycle</b>   |

## Step 5. 2nd Round of Nested PCR

The second PCR amplification step uses the nested GSP2 Anchor Primer to ensure the specific amplification of the correct 3' L1H-specific fragments and add an additional anchor sequence to the ends of the amplified DNA.

1. Dilute the reaction from the First Round of PCR (Step 4 above) 10-fold with water. Use 2  $\mu$ l of the diluted First Round PCR reaction for the Second Round of Nested PCR.
2. Prepare the Second Nested PCR Reaction Mix as follows on ice in a 50  $\mu$ l thin wall PCR tube:

| Component                                 | Volume                      |
|---|-----------------------------|
| 10-fold diluted DNA from 1st PCR (Step 4) | 2 $\mu$ l                   |
| HF DNA Polymerase Buffer (5X)             | 10 $\mu$ l                  |
| dNTP Mix (50X)                            | 1 $\mu$ l                   |
| GSP2 Anchor                               | 1 $\mu$ l                   |
| Adapter Primer                            | 1 $\mu$ l                   |
| PCR-Grade Water                           | 34.5 $\mu$ l                |
| HF DNA Polymerase (100X)                  | 0.5 $\mu$ l                 |
| <b>Total</b>                              | <b>50 <math>\mu</math>l</b> |

3. Mix gently and centrifuge briefly to collect droplets.
4. Amplify the Second Nested PCR Reaction Mix using the following PCR program:

|                  |                  |
|------------------|------------------|
| 98°C, 60 seconds | <b>1 cycle</b>   |
| 98°C, 15 seconds |                  |
| 62°C, 15 seconds | <b>15 cycles</b> |
| 72°C, 30 seconds |                  |
| 72°C, 60 seconds | <b>1 cycle</b>   |
| 4°C, $\infty$    | <b>1 cycle</b>   |

## Step 6. 3rd Round PCR with Indexed Primers

In this amplification step, flanking sequences for the P5 and P7 Illumina adapters are added to the 3' L1 DNA product using the Adapter P5 Primer and one (from the 24) P7 Indexed Primers. Samples with different P7 Indexed Primers can be deconvoluted based on the unique index sequence when loaded in the same lanes of the Illumina Sequencing Cell (see Appendix B: Structure of the Adaptor Sequences of the 3'-Amplicon). As a result, the DNA products with different indexes can be mixed together for cluster generation on the Illumina NGS Sequencing Cell.

1. Use 2  $\mu$ l of the Second Nested PCR product from Step 5 for the final P5/P7 PCR reaction below.
2. For each indexing reaction, use the appropriate Indexing Primer (one of the P7 Ind1-24 primers). Samples that will be sequenced in the same lane of the Illumina cell must have different P7 Index Adaptors. Be certain to keep track of which P7 Index Adaptor was used to amplify each sample so the sequencing results can be deconvoluted to assess the specific reads for each sample.
3. Prepare the reaction mix on ice in a 50  $\mu$ l thin wall PCR tube:

| Component                               | Volume                      |
|---|-----------------------------|
| DNA sample from 2nd PCR (Step 5)        | 2 $\mu$ l                   |
| HF DNA Polymerase Buffer (5X)           | 10 $\mu$ l                  |
| dNTP Mix (50X)                          | 1 $\mu$ l                   |
| P7 Indexed Primer 1-24 (as appropriate) | 1 $\mu$ l                   |
| P5 Primer                               | 1 $\mu$ l                   |
| PCR-Grade Water                         | 34.5 $\mu$ l                |
| HF DNA Polymerase (100X)                | 0.5 $\mu$ l                 |
| <b>Total</b>                            | <b>50 <math>\mu</math>l</b> |

4. Mix gently and centrifuge briefly to collect droplets.
5. Amplify Indexed-3'L1H-DNA products using the following PCR program:

|                  |                |
|------------------|----------------|
| 98°C, 60 seconds | <b>1 cycle</b> |
|------------------|----------------|

|  |                  |
|--|------------------|
| 98°C, 15 seconds<br>62°C, 15 seconds<br>72°C, 30 seconds | <b>10 cycles</b> |
| 72°C, 60 seconds   | <b>1 cycle</b>   |
| 4°C, ∞   | <b>1 cycle</b>   |


6. Analyze the P7-Index-3'L1H-DNA-P5 product on a 3% Agarose TAE gel. Cast a 3% agarose gel with 1X TAE buffer and 400 ng/ml EtBr (60 µg). Mix 5 µl of each PCR product from the Step 6 amplification with 5 µl of 2X loading buffer (e.g. 50% sucrose with Bromophenol Blue tracking dye) and load onto the gel. In at least 1 well, include a 100-bp DNA ladder.
7. Run gel at 5 V/cm until the 100-bp band is about an inch from the bottom. PCR products should appear as a smear of approximately 250-500 bp with the brightest region around 350-450 bp (see Example Gel Data image).

**Optional Analysis Step:** To validate the quality of the Indexed-3'L1H-DNA, use the “Library Validation for Illumina NGS by TA Cloning” procedure described in Appendix A. This may be useful to check if you are having trouble getting good sequencing runs.

## Step 7. Probe Preparation for Illumina NGS

1. Cast a 150 ml 3% agarose gel with 1X TAE buffer and 400 ng/ml EtBr (60  $\mu$ g) using a gel comb that can accommodate approximately 60  $\mu$ l in each well. *Recommended well size:* 1 mm (length) x 8 mm (width) x 7 mm (height).
2. Mix 30  $\mu$ l of each PCR product from the Step 6 amplification with 3  $\mu$ l of 10X loading buffer (e.g. 50% sucrose with Bromophenol Blue tracking dye) and load onto the gel. In at least 1 well, include a 100-bp DNA ladder.
3. Run gel at 5 V/cm until the 100 bp band is about an inch from the bottom (usually about 1 hour), then using a UV Transilluminator and clean scalpel, excise a 2 mm narrow band from the brightest region in the 350-450 bp size range (see Example Gel Data image).

**!** **CAUTION!** Be sure to use UV safety glasses to protect your eyes when viewing and excising DNA from the gel on the UV Transilluminator.



*Example Gel Data*

4. Isolate each PCR product from the gel fragment using the QIAquick Gel Extraction kit and following the manufacturer's protocol. Elute in 30  $\mu$ l of QIAGEN EB Buffer.
5. Measure concentration of each purified P7-Index-3'L1H-DNA-P5 product using a NanoDrop 2000 (or similar spectrophotometer) at OD260. Combine different Index-tagged DNA samples together at



equimolar amounts into the final sequencing sample pool. Dilute the combined DNA samples to 10nM (approximately 2 ng/ $\mu$ l).

## Step 8. Illumina NGS and Data Analysis


---

1. Follow the standard Illumina procedures for Cluster Generation starting with 10 nM of the purified PCR product from Step 6. The final samples can contain up to 24 indexed libraries for sequencing on either the Illumina NextSeq or HiSeq instrument.
2. Add 6  $\mu$ l of each of the custom sequencing primers into the sequencing mix as directed by the Illumina NGS protocol:
  - Read 1 Seq Primer
  - Index 1 Seq Primer
  - Read 2 Seq Primer
3. Perform Illumina NGS Paired-End (PE) sequencing (75 cycles) using the Illumina NextSeq or HiSeq. Proceed to cluster amplification using the Illumina Paired-End Cluster Generation Kit; refer to the manufacturer's instructions for this step. The optimal seeding concentration for cluster amplification of DNA libraries is approximately 6 to 8 pM.
  - Program:

| Step            | Primer             | Cycles |
|-----------------|--------------------|--------|
| <b>Read 1:</b>  | Read 1 Seq Primer  | 60     |
| <b>Index 1:</b> | Index 1 Seq Primer | 6      |
| <b>Index 2:</b> | Read 2 Seq Primer  | 6      |
| <b>Read 2:</b>  | Read 2 Seq Primer  | 13     |

4. NGS FASTQ sequencing files should be analyzed for quality using [FastQC](#). Remove PCR duplicate reads, then map sequences to the appropriate reference genome. During development of this assay, we used reference repetitive element annotation files from the UCSC genome browser and annotated the GRCh38 (hg38) reference genome in CLC Genomics Workbench (CLC Bio/QIAGEN, Cambridge MA). Generate sequencing profiles of L1 insertions from cancer samples and matched normal controls as the reference to isolate all *de novo* cancer-specific L1 insertions.

### Representative Results



*Distribution of L1 elements along the chromosome. The position of each element is covered by several hundred reads on average.*

# Appendices

---

[Appendix A: Library Validation for Illumina NGS by TA Cloning](#)

[Appendix B: Structure of the 3' L1H Amplicon and Primer Locations](#)

[Appendix C: Index Sequences](#)

# Appendix A: Library Validation for Illumina NGS by TA Cloning (Optional Procedure)

This procedure can be used to validate the structure of the amplified fragments and ensure they contain intact P5/P7 sequences for Illumina cluster generation, correct adapter sequences, as well as defined mobile element-chromosome junction regions.

1. Prepare the following reaction mix on ice in a 50  $\mu$ l thin-wall PCR tube:

| Component                            | Volume                      |
|--------------------------------------|-----------------------------|
| DNA sample from 3rd PCR (Step 6)     | 2 $\mu$ l                   |
| Taq DNA Polymerase Buffer (10X)      | 5 $\mu$ l                   |
| dNTP Mix (50X)                       | 1 $\mu$ l                   |
| P7 Indexed Primer 1 (as appropriate) | 1 $\mu$ l                   |
| P5 Primer                            | 1 $\mu$ l                   |
| PCR-Grade Water                      | 39 $\mu$ l                  |
| Taq DNA Polymerase (100X)            | 0.5 $\mu$ l                 |
| <b>Total</b>                         | <b>50 <math>\mu</math>l</b> |

2. Mix gently and centrifuge briefly.
3. Amplify using the following PCR program:

|                  |                |
|------------------|----------------|
| 98°C, 60 seconds | <b>1 cycle</b> |
| 98°C, 15 seconds |                |
| 62°C, 15 seconds |                |
| 72°C, 30 seconds |                |
| 72°C, 60 seconds | <b>1 cycle</b> |
| 4°C, $\infty$    | <b>1 cycle</b> |

4. Follow the instructions in the TA PCR Cloning Kit (Invitrogen/ThermoFisher) to clone 1  $\mu$ l of the PCR reaction into the TA cloning vector and transform the TA cloning mix into the appropriate bacterial strain.
5. Pick 20 colonies at random from the plate. For colony PCR, use the T3/T7 PCR protocol. Sequence the 20 random colonies with T3/T7 primers. If the 3' L1 DNA was prepared correctly, the amplified fragments should contain intact P5 and P7 sequences for Illumina cluster generation and correct adapter sequences, as well as defined mobile element-chromosome junction regions as shown on Figure 2.

# Appendix B: Structure of the 3' L1H Amplicon and Primer Locations


## A. Consensus Sequence of L1H 3' UTR

5' - GGAGGGATAGCATTGGGAGATATACCTAATGCTAGATGACACATTAGTGGGTGCAGCGCACCAGCATGGCACATGTATACATATGTA  
GGAGGGATAGCATTGGGAGAT -----> 3'L1H GSP1      3'L1H GSP2 CGCACCAGCATGGCACATGT ----->

ACTAACCTGCACAATGTGCACATGTACCCTAAAACCTTAGAGTATAATAAAAAAAAA(A)<sub>n</sub>---Insertion site---Chromosome DNA-3'

(A) Consensus sequence of the L1H 3' UTR portion amplified in the protocol.

## B. 3' L1H Amplicon Structure



(B) 3'L1H amplicon structure.

**Structure of 3' L1H Amplicon. (A) Consensus sequence of the L1H 3'UTR portion amplified in the protocol.** Primers specific for L1 sequences are used for amplification of only the youngest families of L1 (L1H) elements from the human genome. Two primers (Gene-Specific Primers, or GSP1 and GSP2) are designed for 3' ends of the elements. For L1H, 3' L1H UTR-specific primers bind within the UTRs of the L1 and extend away from the element and into the chromosome insertion sites. The primers are designed in such a way that they specifically amplify only the youngest subfamilies of polymorphic transposable elements which are still active in the genome, as opposed to the ancient elements fixed in the population. **(B) 3' L1H amplicon structure.** 3' L1H fragments are amplified by 1st Round PCR followed by 2nd Round Nested PCR in Step 5. Illumina NGS P5 and P7 adapter sequences (green) are introduced during the 3rd PCR in Step 6. The library is then sequenced by Illumina NGS in Step 8 using the Read 1 Seq Primer to read 13 nt of the L1H, and 3' L1H Read 2 Seq Primer to read the corresponding chromosome junction portion as well as the Index 1 Seq Primer to read the multiplexing index shown in red.

# Appendix C: Index Sequences

---

## Index Sequences for P7 Indexed Primers 1-24

| Index | Sequence |
|-------|----------|
| 1     | CGTGAT   |
| 2     | ACATCG   |
| 3     | GCCTAA   |
| 4     | TGGTCA   |
| 5     | CACTGT   |
| 6     | ATTGGC   |
| 7     | GATCTG   |
| 8     | TCAAGT   |
| 9     | CTGATC   |
| 10    | AAGCTA   |
| 11    | GTAGCC   |
| 12    | TACAAG   |
| 13    | TTGACT   |
| 14    | GGAACT   |
| 15    | TGACAT   |
| 16    | GGACGG   |
| 17    | CTCTAC   |
| 18    | GCGGAC   |
| 19    | TTTCAC   |
| 20    | GGCCAC   |
| 21    | CGAAAC   |
| 22    | CGTACG   |
| 23    | CCACTC   |
| 24    | GCTACC   |



# Technical Support

---

Please contact us by phone or email if you have any questions.

Phone: +1 (650) 938-3910

Toll-Free: +1 (877) 938-3910

Fax: +1 (650) 938-3911

**E-mail:**

Technical Support: [tech@cellecta.com](mailto:tech@cellecta.com)

General Information: [info@cellecta.com](mailto:info@cellecta.com)

Sales: [sales@cellecta.com](mailto:sales@cellecta.com)

Orders: [orders@cellecta.com](mailto:orders@cellecta.com)

**Postal Mail:**

Cellecta, Inc.

320 Logue Ave.

Mountain View, CA 94043

For the latest technical news and updates, visit Cellecta's blog at: <https://www.cellecta.com/blog-news/>

# Contact Us

---

## Email Addresses

General information: [info@cellecta.com](mailto:info@cellecta.com)

Sales: [sales@cellecta.com](mailto:sales@cellecta.com)

Orders: [orders@cellecta.com](mailto:orders@cellecta.com)

Technical Support: [tech@cellecta.com](mailto:tech@cellecta.com)

## Mailing Address

Cellecta, Inc.  
320 Logue Ave.  
Mountain View, CA 94043  
USA

## Phone Numbers

Phone: +1 650 938-3910

Toll-free: (877) 938-3910

Fax: +1 650 938-3911


For the latest technical news and updates, visit Cellecta's blog at:

<https://www.cellecta.com/blog-news/>



# PDF Download Link

---

 [Download the PDF Version of the LINE-1 Mobile Element NGS Assay User Manual here](#)