



Design of shRNA libraries in pMSCV vectors with the HTS6M cassette.

The GG214-GT-C7 HTS6M Library Cassette Design allows standard two-round amplification of the shRNA cassette from genomic DNA. In the first step, the cassette is amplified with F2 + R2 primers from genomic DNA. In the second step, the enriched cassette is amplified with Gex1-NF2 + Gex2-NR2 primers (Standard Next-Gen Sequencing) and PCR products are sequenced on the HiSeq 2000 or 2500. Amplified PCR products (approx. 200bp) are compatible only with single read (SR) Illumina flow cells.

Cassette design: 26M6-HTS6M
Oligo size: 139n
shRNA Design: GG214-GT-C7
NGS cassette: Version 6M

Related Vectors: pMSCVURP
NGS Platforms: HiSeq 2000/2500, GAIIX
(Please contact Cellecta at tech@cellecta.com for NextSeq primers)

Visit the Cellecta website for User Manuals and Vector information:
User Manuals, PACs: <https://www.cellecta.com/resources/product-manuals-and-certificates/>
Cellecta Vectors: <https://www.cellecta.com/resources/vector-information/>

BpiI Sense (21) -GT Stem4-LoopC7-Stem4 Antisense (21) EcoRI shRNA Barcode (22)
TCA GAAGAC GCACCGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNTTAAATATTCATAGCNNNNNNNNNNNNNNNNNNNNNNNTTTTGAATTCGCACCAGCAGCTACGCANNNNNNNNNNNNNNNNNNNNNNTTCGCC GICTTCGT
BpiI

shRNA Construct/Library in pMSCV vector:

ClaI-FwdU6-1>FwdU6-3>U6-shRNA-T6-F2>Gex1-NF2>shRNA BC22<GexSeqS-cPPT<Gex2-NR2<R2<RevUbiC1

F2

Gex1-NF2 5' -TCGGATTTCGACCAGCAGCTA

5' -TCAAGCAGAAGACGGCATAACGATCGCACCAGCAGCTACGCA XhoI Stuffer

U6 -ACCGG-shRNA-TTTTGAATTCGCACCAGCAGCTACGCA-BC22-TTCGGACTGTAGAAGTCTGAACTCAGCAATAAAAGAATGGTACAGTGCAGAAACAATAGTAGACATAATAGCAACAGA

U6 -TGGCC-shRNA-AAAAACTTAAGCGTGGTCGTGCGATGCGT-BC22-AAGCCTGACATCTTGAGACTTGAGCTCGTTATTTTCTTACCATGTACAGTCTTTGTTATCATCTGTATTATCGTTGTCT

AAGCCTGACATCTTGAGACTTGAGCA-5' GexSeqS

EcoRI AgeI BstBI UbiC promoter

CATACAACTAAAGAACTGCGTTGTTGTCGGTGCCTGTTCTCTGCTCTTCACGCTACTGAATTCATCACCGGTTTCGAAAGGCCTCCGCGCCGGGTTTTGGCGCCTCCCGGGGGCCTTCCTCAC

GTATGTTTGAATTCCTTAGACGCAACAACAGCCAGGACGAGCAAGAGACGAGAAGTGCATGACTTAAGTAGTGGCCAAGAAGCTTCGGGAGGCGCGCCCAAACCGCGGAGGCGCCCGGGGGGAGGAGTG

AGACGCAACAACAGCCACGAGAGCCACCAGCGGCATAGTAA-5' AGGAGTGCCGCTCAGCGCGGA RevUbiC1

Gex2-NR2 AAGAGACGAGAAGTGCATGA-5' R2



Amplification of shRNA-specific Barcodes using two rounds of amplification (starting from 100-200ug of genomic DNA with Clontech's Titanium Taq DNA polymerase)

First Round

F2 TCGGATTCGCACCAGCAGCAGCTA
R2 AGTAGCGTGAAGAGCAGAGAA

Second Round - amplicon size ~200bp

Gex1-NF2 TCAAGCAGAAGACGGCATAACGATCGCACCAGCAGCTACGCA
Gex2-NR2 AATGATACGGCGACCACCGAGAGCACCAGACAACAACGCAGA

Sequencing Primer (22 cycles)

GexSeqS AGAGGTTCAGAGTTCTACAGTCCGAA

pMSCV Standard NGS Amplification/Sequencing Primers:

| Primer Name | Used for | Sequence (IDT preferred) |
|-------------|-----------------------|--|
| F2 | 1 st Round | 5' - TCGGATTCGCACCAGCAGCAGCTA - 3' |
| R2 | 1 st Round | 5' - AGTAGCGTGAAGAGCAGAGAA - 3' |
| Gex1-NF2 | 2 nd Round | 5' - TCAAGCAGAAGACGGCATAACGATCGCACCAGCAGCTACGCA - 3' |
| Gex2-NR2 | 2 nd Round | 5' - AATGATACGGCGACCACCGAGAGCACCAGACAACAACGCAGA - 3' |
| GexSeqS | NGS | 5' - AGAGGTTCAGAGTTCTACAGTCCGAA - 3' (HPLC Purified) |

Standard Insert Screening PCR Primers:

FwdU6-1 CAAGGCTGTTAGAGAGATAAATTGGAA
Rev-cPPT10 TGTATGTCTGTTGCTATTATGTCTAC

Standard Sanger Sequencing Primer

FwdU6-3 ATTAGTACAAAATACGTGACGTAGAA